

Kamoer Fluid Tech (Shanghai) Co., Ltd.

Version: A/2

Product name	Adjustable peristaltic pump		
Product model KCP200			
Execution date	2023.09.22		
Company name	Kamoer Fluid Tech (Shanghai) Co.,Ltd.		



Model: KCP200

File No .: CPBZ-KCP200-01

Version No.: A2


Effective Date: 2023-09-22

1. Product Overview

1.1 Product physical picture



1.2 Product assembly drawing



1.3 Model definition (stepper motor drive)

	KCP200-B16-GB	KCP200-B16-US	KCP200-B16-UK	KCP200-B16-EU	KCP200-B16-AU
Model	National standard	American	British Standard	European standard	Australian
		Standard			Standard
Sales code	CK.33.18.0001	CK.33.18.0002	CK.33.18.0003	CK.33.18.0004	CK.33.18.0005

Add Building , No. XXiangjing Road, Songjiang District, Shanghai







Model: KCP200

File No .: CPBZ-KCP200-01

Version No.: A2

Effective Date: 2023-09-22

1.4 Pump tube material

Pump tube material	Characteristic	Standards compliant
Silicone tube (S)	Low adsorption, good temperature resistance, good chemical corrosion resistance, life expectancy of 200 hours	RoHS FDA
BPT tube (B)	Resistant to common chemicals, acid and alkali, and oxidation resistance, with a lifespan of 1000 hours	RoHS

1.5 Performance characteristics and typical applications

1.5.1 Performance characteristics

- Low-cost adjustment of flow rate
- Easy and quick replacement of pump tubing
- Low noise, small space occupation
- Simple structure, maintenance-free
- With three rotors, moderate pulsation

1.5.2 Typical application

- Laboratory: liquid transfer, reagent packing
- Industrial process: such as pumping hydraulic oil
- Fine chemicals: such as cosmetic packaging, slurry recovery and filtration
- Food Process: Beverage Filling







Ion Chromatograph

Auto Car Wash

Five classification blood analyzer

Hypochlorite water generator

Add Building , No. Wiangjing Road, Songjiang District, Shanghai



File No .: CPBZ-KCP200-01

Version No.: A2

Name: peristaltic pump | Model: KCP200 | Effective Date: 2023-09-22

1.6 Product main material

Chassis	KPHM200 pump head			
Cilassis	Cover	Sync disk	Rotor	Fitting
ABS plastic	PA6	РОМ	PET	PP

1.7 Product Certification, Intellectual Property

CN 307203632S




1.8 Product Risk Warning

- Hoses are resistant to liquid media Need to check chemical compatibility or conduct immersion tests, improper hose selection may cause hose damage quickly
- The working environment of the product should not exceed 45°C, and the humidity should not exceed 70% (no condensation). The harsh working environment will cause premature damage to the product
- Fluid leakage accidents due to hose ruptures depend on the fluid medium and your specific application conditions
 - High overloading may result in premature failure of the product

2. Product specification

2.1 Product Size



Add Building , No. XXiangjing Road, Songjiang District, Shanghai

Tel ⊠ -XXXX 🛛





Model: KCP200

File No .: CPBZ-KCP200-01

Version No.: A2

Effective Date: 2023-09-22

2.2 Technical Parameters

	Project	Technical Parameters			
	Using electric	DC24V 1 A power adapter			
	Power	<24W			
	Drive	DC brushed motor			
	Working	Toron outstand 0 - 40°C valeting bound dita < 000/ DII			
	environment	Temperature 0∼40°C, relative humidity<80%RH			
	Total Weight	About 0.52 Kg			
Ва	Reference noise	\leq 60dB (test environment noise 34dB , the horizontal distance between the test product and			
sic p	value	the noise meter is 0. 5 m)			
Basic parameters	Range of rotation	<400rpm			
mete	Speed adjustment				
SLS	resolution	/			
	Flow range	100~200ml/min			
	Display method	/			
	Control method	Potentiometer			
	Operating mode	Speed regulation mode			
	Way of working	manual			
	Support function	Start and stop, speed regulation			
		Technical Parameters			
	Model code	B16 (16#)			
Pu	Pump tube	DI LOCT			
Pump flow	material	Pharmed BPT			
low	ID x OD	Metric: 3.2 mm x 6.4 mm Inch: $1/8$ in $\times 1/4$ in			
	Number of rotors	3			
	Max flow	200ml / min			
Precautions		The maximum flow rate above is obtained at room temperature (about 25 °C), using a new pump tube that has been aged for 30 minutes and tested water at 400 rpm, for reference only. Ambient temperature, the material and elasticity of the pump tube, the viscosity of the test liquid and other factors will affect the actual flow rate. The thickness of the pump tubing will affect the actual maximum speed for stable operation.			





File No .: CPBZ-KCP200-01

Version No.: A2

Effective Date: 2023-09-22

2.3 Accessories

Assembly name	Accessories picture		Function introduction
Connect pump tubing 14.03.0036		3.2mm × 6.4mm	Connect the transmission between the liquid source and the pump, low adsorption, good temperature resistance, good chemical corrosion resistance, passed RoHS FDA.

